

**Protocol for fluorescence-based neuraminidase assay and
neuraminidase inhibition assay
(St. Jude Children's Research Hospital)**

Buffers and solutions:

0.1 M MES (pH=6.5)

2.1 g MES (*Sigma, M5287-250G*)

100 ml H₂O

Check pH (adjust with conc. NaOH)

MES hydrate > 99.5% titration, C₆H₁₃NO₄SxH₂O,

FW 195.2 (Anh)

325 mM MES (pH=6.5)

6.9 g MES

100ml H₂O

Check pH (adjust with conc. NaOH)

0.1 M CaCl₂

1.47 g CaCl₂x2H₂O

100 ml H₂O

10 mM CaCl₂

0.15g CaCl₂x2H₂O

100 ml H₂O

Enzyme buffer

32.5 ml 0.1 M MES (pH=6.5)

4.0 ml 0.1 M CaCl₂

63.5 ml H₂O

Store at +4°C

Standard

4 mM Stock standard solution (stable for ~ 1 month at 4°C, protect from light)

7 mg standard (4-Methyl-Umbelliferone) (*Sigma; Cat= M-1381*)

5 ml absolute ethanol

5 ml 0.9 % NaCl (wt/v)

40 μM Working standard solution (prepare right before experiments)

10 μl stock standard solution

990 μl 0.9 % NaCl (wt/v)

Stop solution (pH=10.7)

100 ml ethanol

250 ml H₂O

Mixed water and ethanol first and then add Glycin

3.004 g Glycin (*Fisher; G46-500*)

Check pH (adjust with conc. NaOH)

Add water to 400 mls

Preparation of substrate stock (5 mM solution):

Substrate: 2'-(4-Methylumbelliferyi)- α -D-N-acetylneuraminic acid, sodium salt
Sigma-Aldrich M-8639, MW 489.4

1. Prepare water-saturated methyl-tert butyl ether (*Methyl-tert butyl ether, HPLC grade; e.g., Fisher E127*) by shaking 20 ml ether with 10 ml distilled H₂O in glass-stoppered glass Erlenmeyer (KILMAX #22, 125 ml volume, Fisher Scientific). Remove (lower) water phase after separation of water and ether layers and shake again with another 10 ml dH₂O. Remove this water as well. All procedures must be done with glass pipets or pateur glass pipets.
2. Add 2.25 ml distilled H₂O to a 5 mg vial of substrate [2'-(4-Methylumbelliferyi)- α -D-N-acetylneuraminic acid, sodium salt]. It is better to order the substrate that is aluquited by 5 mg per vial/bottle; in this case it is possible to add distilled H₂O directly to the vial/bottle. The powder dissolves quickly. Transfer solution to a 15 ml polypropylene snap-cap tube or a glass tube.
3. Add 5 ml water-saturated methyl-tert butyl ether to substrate solution. Shake vigorously. Remove (upper) ether phase by paster pipet. Repeat this twice. In total, three extractions.
4. Put substrate solution on ice. To remove the remaining ether, nitrogen gas is bubbled through the substrate solution for ± 2 hours, until the solution has lost almost all ether smell.
5. Aliquot substrate into 100 or 50 μ l per tube and store at -20°C.

Substrate mixture for NA assay

Prepare the mixture just before the experiment and leave on ice

Concentration of substrate mixture – 167 μM , after adding to the plate – the final concentration of substrate will be 100 μM .

Buffers/Substrate	½ plate	1 plate	2 plate	3 plate
1 Volume 325 mM MES	250 μl	500 μl	1000 μl	1500 μl
3 volume 10 mM CaCl_2	750 μl	1500 μl	3000 μl	4500 μl
Substrate	50 μl of 5 mM MUN (from -20C) 450 μl H ₂ O	100 μl of 5 mM MUN (from -20C) 900 μl H ₂ O	200 μl of 5 mM MUN (from -20C) 1800 μl H ₂ O	300 μl of 5 mM MUN (from -20C) 2700 μl H ₂ O

Preparation of the stocks of NA inhibitors (Compounds to be used in NI assays)

Molecular weight:

Oseltamivir carboxylate, tartrate salt: 386.4 g/mol

Zanamivir monohydrate: 350.3 g/mol

Note: It is important to request the compounds from an appropriate source and ask about Product information. The composition and purity varies.

Steps:

1. In our experiments we are using the stock of NA inhibitor at concentration **2mM**.

For oseltamivir carboxylate (active ingredient of oseltamivir):

- Weight 0.01544 gram oseltamivir carboxylate and dissolve in 20 mL of sterile ddH₂O.

- Explanations:

1M (1000mM) – 386 gram – 1L (1000 mL)

2 mM - X gram – 1L (1000 mL)

X= 0.772 gram per 1L (1000 mL) or 0.015444 gram per 20 mL

For zanamivir:

- Weight 0.014 gram zanamivir and dissolve in 20 mL of sterile ddH₂O.

- Explanations:

1M (1000mM) – 350 gram – 1L (1000 mL)

2 mM - X gram – 1L (1000 mL)

X= 0.70 gram per 1L (1000 mL) or 0.014 gram per 20 mL

2. Vortex the mixture and the compound can be easily dissolved. The optional procedure is to incubate the compound in water bath at 37C for 15-20 minutes.
3. Filter the solution through a sterile Acrodisc Syringe Filter (0.2 µm HT Tuffryn Membrane Low Protein Binding). This procedure is performed in a hood with a syringe and allows elimination of bacteria and fungus.
4. Aliquot the NA inhibitor solution (at a concentration 2 mM) into 200µL per tube in sterile 0.5 mL tubes. This procedure is performed in a hood.
5. Store at -20°C freezer and avoid freeze-thawing.

NA inhibition assay

- Mix 10 μl virus in enzyme buffer (1500 to 2000 relative fluorescent units, as determined in NA assay) with 10 μl compound in enzyme buffer. Different concentrations of the compound must be used. Incubate the plate at 37° C for 30 minutes (with shaking sometimes).

Design of the plate

Blank	Standard 5 μl = 200pmol	0.5 μM + Virus 1	0.1 μM	0.05 μM	0.01 μM	0.005 μM	0.001 μM	0.0005 μM	0.0001 μM	0.00005 μM	virus + enzyme buffer
Blank	Standard 5 μl = 200pmol	Virus 1	→	→	→	→	→	→	→	→	→
Blank	Standard 10 μl = 400pmol	Virus 2									
Blank	Standard 5 μl = 400pmol	Virus 2									
Blank	Standard 20 μl = 800pmol	Virus 3									
Blank	Standard 20 μl = 800pmol	Virus 3									
Blank	Empty	Virus 4									
Blank	Empty	Virus 4									

- Add 30 μl per well of the substrate mix (final concentration of the substrate 100 μM) and further incubate 1 hour at 37° C.
- The reaction is terminated by the addition of 150 μl per well of stop solution (Glycin in ethanol, pH=10.7).
- The fluorescence of the released 4-methylumbelliferone was measured in a Fluoroskan II, type 371 (Labsystems, Helsinki, Finland) spectrophotometer using excitation and emission wavelengths of 355 and 460 nm, respectively.

REFERENCE:

Potier M, Mameli L, Belislem M, Dallaire L, Melanxon SB (1979) Fluorometric assay of neuraminidase with a sodium (4-methylumbelliferyl- α -D-N-acetylneuraminate) substrate. Anal Biochem 94: 287-296.